

Environmental control Documentation for **EFP 2005 doorphone**

To be used together with Picolo, PicoCombi and Erifox

Revision history

When changes are made on the product or in the documentation, this will be given a revision number. The list below keeps track of these changes.

Date	Name	Description	
19.10.05	Bent-Håvard Sollid.	First version with English text.	
01.07.06	Bent-Håvard Sollid.	EC Declaration of Conformity updated with new	
		directives.	
		Added WEEE icon.	
		Changed revisionlist to new type.	
01.08.08	Bent-Håvard Sollid.	New version of outdoor unit:	
		- Ch 6.2.1, changed text.	
		- Ch, 6.2.3, 6.2.4, changed wiring diagram.	
01.01.12	Bent-Håvard Sollid.	Added wiring diagram for the new version of PicoCombi,	
		MK II.	

Contents – Documentation EFP 2005

1	Intr	oduction	2
2	Get	ting started, installation	2
3		r manual	
	3.1	Accessories	
4	Mai	ntenance	
	4.1	Cleaning	
	4.2	Outdoor use	
5	Tro	ubleshootingubleshooting	3
	5.1	If errors arise	
6	Tecl	hnical information	
	6.1	Fixed mounting	
	6.2	Electrical connections	4
	6.2.	1 Outdoor unit	4
	6.2.2	2 Amplifier box	4
	6.2.3	Wiring diagram for Picolo/PicoCombi version I	6
	6.2.4		
	6.2.5		
	6.3	Default (recycling)	
	6.4	Technical data	
	6.5	CE markings	



1 Introduction

This documentation is made for Picomeds doorphone EFP 2005, which is a part of Picomeds environmental system. It consists of 2 units:

Amplifier-box made for indoor placement. It consists of necessary amplification and voicedirection control.

Outdoor unit consisting of a box with a speaker, a microphone and a push-button.

From the amplifier box there shall be cables to:

Indoor unit (Picolo, PicoCombi or Erifox PABX). Signals and power.

Outdoor unit. Signals.

The equipment does not have any handset, only loudspeaking function. There is an built-in voice controlling function which controls voicedirection, i.e if the loudspeaker or the microphone shall be enabled. It is possible to regulate the volume from the outdoor unit, but it is recommended that the preset value is kept unaltered. EFP 2005 shall be powered with 12 VDC.

2 Getting started, installation

The following parts are included when the product is delivered:

Amplifier box, made to be connected to Picolo/PicoCombi or Erifox. A Picolo/PicoCombi version can not be connected to an Erifox version or vice verca.

Outdoor unit.

This documentation.

The equipment is developed to be connected to Picomeds loudspeaking indoor units Picolo and PicoCombi, or an Erifox PABX. Indoor unit is not included in the EFP 2005 shipment.

EFP 2005 made for Picolo/PicoCombi can not be connected to an Erifox version or v.v. The Erifox version is clearly marked with "Erifox".

Find a suitable place to mount the outdoor unit. It is recommended to place it where it is shielded for wind and precipitation. There is space for the users name on the outdoor unit. It can be marked with a label, waterproof pencil etc.

Find a suitable place to mount the amplifier box. It shall be mounted indoor in a dry place. It is recommended that it is not more than 4-5 meters cable between the outdoor and the amplifier box.

Mount cable between amplifier box and outdoor unit, connect as shown in the wiring diagram.

Mount cable between amplifier box and indoor unit, connect as shown in the wiring diagram. Erifox version: the telephone line used for doorphone has to be enabled. This is done in programming mode. Find details in the PABX usermanual.

The doorphone is now ready to be used together with an indoor unit. Test the communication. The outdoor units volume in the speaker can be adjusted with the potentiometer marked with "GAIN" on the PCB inside the amplifier box. It is recommended that the preset value is kept unaltered. A red light indicator (LED) is illuminated when the amplifier box is powered with 12 VDC. Position for the "GAIN" potentiometer and the LED is shown in chapter 6.2.3.

3 User manual



The outdoor unit has a pushbutton which shall be pressed by the visitor. The conversation is enabled by the user from the indoor unit (Picolo, PicoCombi or telephone via PABX). This process is described in the indoor units documentation.

3.1 Accessories

Accessories for the EFP 2005 doorphone.

The outdoor units frontfoil can be replaced if the users name shall be changed.

4 Maintenance

EFP 2005 doorphone is developed to work without any need for adjustment etc. There are no fuses or batteries which has to be changed.

4.1 Cleaning

The equipment should be cleaned when necessary. Wipe the equipment with a slightly damp cloth or an antistatic cloth. Never use a dry cloth, while this result in risk of static charges. It shall not be exposed for fluids in any way.

4.2 Outdoor use

Picomeds EFP 2005 doorphone shall have the amplifier box mounted indoor, while the outdoor unit should be mounted in a shielded place. Attempt to find a place where it not will be exposed by wind, rain and snow. The reason for this, is that water will be able to get inside the box and degrade the equipments function.

5 Troubleshooting

If your equipment does not behave as expected, first of all try to solve the problem yourself using the following table:

Symptom	Possible cause	Recommended action
Completely dead. Light	No power connected to the	Ensure that the power supply
indicator not illuminated.	indoor unit (the EFP 2005 is	is the correct one which was
	powered from the indoor	delivered with the equipment.
	unit).	Ensure that the power supply
		is connected to a functioning wall outlet.
		Ensure that the power supply
		is properly connected to the
		equipment.
No sound in the speaker in the	<u> </u>	Check the "LS" cable.
outdoor unit.	and outdoor unit.	
No sound in the speaker in the	<u> </u>	Check the "MIC" cable.
indoor unit.	and outdoor unit.	
No sound neither in the	Cable between the amplifier	Check wiring: LS, Mic, La
outdoor unit or in the indoor	box and outdoor unit.	and Lb.
unit.	Cable between the amplifier	
	box and indoor unit.	
Push-button does not	Cable between the amplifier	Check wiring both places.
function.	box and outdoor unit.	
	Cable between the amplifier	
	box and indoor unit.	



There are no fuses inside the device which can be replaced.

5.1 If errors arise

If there are errors which are not possible to solve, contact your supplier of the equipment, or the national distributor. Please make a note beneath of the telephone number to your supplier.

Name	Telephone number, address etc.
Supplier:	
National distributor:	

6 Technical information

The amplifier box is made by:

Injection moulded box in recyclable ABS-94HB.

Printed circuit board (PCB) screwed into bottom part. The PCB has 3 terminal blocks: one for connection to outdoor unit, one for connection to indoor unit and one marked with "FREE".

This is made as a help to the technician if there is a need for joining wires etc. No one of the connecting points are connected to each other or to components on the PCB.

A red LED illuminates when power is connected.

Potentiometer ("GAIN"). Can be used to regulate speaker volume in outdoor unit.

Outdoor unit:

Front part with microphone, loudspeaker and pushbutton.

Frontfoil glued on front part. Picomeds logo, product name and space made for the users name on it with a label, waterproof pen etc.

Backside has 6 electrical connections.

Frame for wall mounting.

6.1 Fixed mounting

The amplifier box can be mounted vertically or horizontally as needed.

The outdoor unit is made to placed vertically and shielded for wind, rain and snow.

6.2 Electrical connections

6.2.1 Outdoor unit

On the back there is terminal blocks with connections. Connect as advised on wiring diagram. Connections are made for maximum 2,5 mm² cable.

Remove 6 mm of the insulation.

6.2.2 Amplifier box

There are 3 terminal blocks in the amplifier box and each of them have 6 connections.

Connections are made for maximum 2,5 mm² cable.

Remove 6 mm of the insulation.

Terminal block at the lower right is for cable to indoor unit.



Terminal block at the lower left is for cable to outdoor unit.

Terminal block at the top left is marked "FREE". This is made as a help to the technician if there is a need for joining wires etc. No one of the connection points are connected to each other or to components on the PCB.

Description of the connections:

12V and GND is powersupply to the amplifier. Only use direct current (DC).

La, Lb is sound between indoor unit and amplifier box.


DB is DoorBell, one-pole switch in outdoor unit. Connects GND back to indoor unit.

LS is LoudSpeaker, sound to loudspeaker in outdoor unit. Non-polarized.

Mic+ and Mic- is sound from the outdoor units microphone. Polarized.




6.2.3 Wiring diagram for Picolo/PicoCombi version I






6.2.4 Wiring diagram for Picolo/PicoCombi version II

OBS: this is for PicoCombi MK II – not same diagram as for MKI.





6.2.5 Wiring diagram for Erifox PABX





6.3 Default (recycling)

The following shall be done on Picomeds EFP 2005 when it shall be used by a new user. If necessary, make a copy of this page, fill inn serial number, version and tick off the points which are done with signature and date at each point, and file it as your routines requires. Serial number:

Check the equipment. Is there any physical damages, cracks etc in box or foil?

Check that foil is not damaged, that it is fastened all around the edges and that it not has loosened from its surface. If it is marked with former users name, remove this marking. If not possible to remove name, replace foil (can be ordered from supplier).

Clean the device.

Ensure that all components are available, amplifier box and outdoor unit.

Make a function check.

The equipment has to be connected to an indoor unit. It is recommended to do this before the equipment is disconnected.

Ensure that all functions are working properly: check if "ding-dong" is responding, and that it is possible to have 2-way communication between out- and in indoor units.

Adjust potentiometer for outdoor loudspeaker volume ("GAIN") clockwise until end-stop, thereafter counter clockwise 90°.

Store.

6.4 Technical data

Type Picomed EFP 2005 doorphone. Operating Pushbutton on outdoor unit.

The conversation is activated from indoor unit (not included).

Voice direction is electronically controlled.

Power supply From indoor unit, 12 VDC.

Current consumption < 250 mA.

Temperature Amplifier box: 0° C to $+40^{\circ}$ C.

Outdoor unit: -20°C to +40° C.

Measurements Amplifier box: LxWxH: 117 x 94 x 14 mm.

Outdoor unit: LxWxH: 130 x 100 x 45 mm.

Weight Amplifier box: 160 g.

Outdoor unit: 250 g.

Programmable parameters Outdoor volume can be adjusted by technician. It is recommended

that the preset value is kept unaltered.

Waste management Dispose as electronic waste.

Developed and manufactured in Norway.

6.5 CE markings

Picomeds environmental control systems are developed and manufactured in accordance with EU directives, national regulations and harmonised standards which are applicable, and thereafter marked with the CE-mark. The Declaration of Conformity is in English version and it is set out on the next page.









EC DECLARATION OF CONFORMITY

Manufacturers name and address:

Picomed as.
Brokelandsheia.
N-4993 SUNDEBRU, Norway.

Declare under our sole responsibility that the product:

EFP 2005

Which is a doorphone consisting of an amplifier (indoor unit) and an outdoor unit.

Compiles with the requirements of the following directives:

- Low Voltage Directive 73/23/EEC.
- EMC Directive 89/336/EEC and 93/68/EØF.
- Medical Devices Directive 93/42/EEC (Class I).
- RoHS Directive 2002/95/EC.
- WEEE Directive 2002/96/EC.

It is in conformity with the following standards:

- EN 60601-1:1988-12 (IEC 601-1:1988):
 Medical electrical equipment Pt. 1: General requirements for safety. Defined as a Class I product.
- EN 60601-1-2 (EMC).
- EN 61000-3-2 (radiation), EN61000-3-3 (radiation) and EN61000-4-3 (Immunity).
- NS-EN 12182:1999: Technical aids for disabled person. General requirements and test methods.
- NS-EN ISO14971: 2001-04: Medical devices. Application of risk management to medical devices (ISO 14971:2000)(Corrigendum AC:2001 incorporated).
- NS-EN ISO 14971:2000/A1: 2003-05: Amendment A1. Medical devices. Application of risk management to medical devices. Amendment 1: Rationale for requirements (ISO 14971:2000/AM1:2003).

Supplementary information:

The products carries the CE mark accordingly. The CE mark was first applied in the year 2005. The technical documentation relevant to the above equipment will be held by us.

Gjerstad, Norway

Location

01. july 2006.

Issue Date

Bent-Håvard Sollid.

Development Manager